

# LATERAL FLOW TESTS

---

POWERFUL IMMUNO DIAGNOSTIC  
TOOLS IN YOUR HANDS



[www.operon.es](http://www.operon.es)



## PRODUCT LIST

### GASTROINTESTINAL

ANALYTE	PRODUCT NAME	SAMPLE
Rotavirus	Simple/Stick Rotavirus	Stool
Rotavirus-Adenovirus	Simple/Stick Rota-Adeno	Stool
Noravirus (GI, GII)	Simple/Stick Noravirus	Stool
Rotavirus-Noravirus	Simple/Stick Rota-Noro	Stool
Rotavirus-Adenovirus-Noravirus (GI, GII)	Simple Rota-Adeno-Noro	Stool
<i>G. lamblia</i>	Simple/Stick Giardia	Stool
<i>C. parvum</i>	Simple/Stick Crypto	Stool
<i>G. lamblia-C. parvum</i>	Simple/Stick Crypto-Giardia	Stool
<i>E. histolytica</i>	Simple/Stick Entamoeba	Stool
<i>G. lamblia-C. parvum-E. histolytica</i>	Simple/Stick Crypto-Giardia-Entamoeba	Stool
<i>C. difficile</i> , GDH antigen	Simple/Stick GDH	Stool
<i>C. difficile</i> , toxins A and B	Simple/Stick 2A-Bdiff	Stool
<i>C. difficile</i> , GDH antigen and toxins A/B	Simple GDH-Toxins	Stool
<i>E. coli</i> verotoxin	Simple/Stick VER-0157	Stool
<i>H. pylori</i>	Simple/Stick <i>H. pylori</i> *	Stool
Calprotectin	Simple Calprotectina	Stool

\*Only available for certain countries

### RESPIRATORY

ANALYTE	PRODUCT NAME	SAMPLE
Respiratory Syncytial Virus (RSV)	Simple RSV	Nasal/Nasopharyngeal swab Nasopharyngeal wash/aspirate
Influenza	Simple/Stick Flu A+B	Nasal/Nasopharyngeal swab Nasopharyngeal wash/aspirate
Adenovirus	Simple/Stick RespirAdeno	Nasal/Nasopharyngeal swab Nasopharyngeal wash/aspirate
RSV-Influenza	Simple RSV - Flu A+B	Nasal/Nasopharyngeal swab Nasopharyngeal wash/aspirate
RSV-Adenovirus	Simple RSV - RespirAdeno	Nasal/Nasopharyngeal swab Nasopharyngeal wash/aspirate
<i>L. pneumophila</i>	Simple/Stick Legio Pneumo	Urine
<i>S. pneumoniae</i>	Simple/Stick Strep Pneumo	Urine
<i>S. pneumoniae-L. pneumophila</i>	Simple/Stick Strep Pneumo-Legio	Urine
Streptococcus, group A	Simple/Stick Strep A	Throat swab

# PRODUCT LIST

## TROPICAL

ANALYTE	PRODUCT NAME	SAMPLE
T. cruzi	Simple/Stick Chagas	Serum
T. cruzi	Simple/Stick Chagas WB	Blood

## FOOD INTOLERANCES

ANALYTE	PRODUCT NAME	SAMPLE
Celiac disease (tTG)	Stick CD1	Serum
Celiac disease (tTG and gliadins)	Simple CD1 WB	Blood
Celiac disease (tTG and gliadins)	Simple CD2 WB	Blood

## TUMOR MARKERS

ANALYTE	PRODUCT NAME	SAMPLE
Prostate specific antigen	Simple/Stick PSA	Serum

## FOOD ANALYSIS

ANALYTE	PRODUCT NAME	SAMPLE
Milk adulteration	Stick cGMP	Milk

## OTHERS

ANALYTE	PRODUCT NAME	SAMPLE
Pregnancy	Simple/Stick hCG s&u	Serum-urine
Hepatitis B	Simple/Stick HBsAg*	Serum/Plasma

\*Only for sales in non-EU countries

# Infectious diseases *Enfermedades infecciosas*

GDH, 2 A-B diff

Immunochromatographic test for the diagnosis of Infectious diseases in stool samples  
Test inmunocromatográfico para el diagnóstico de enfermedades infecciosas en heces



**operon**  
immuno & molecular  
diagnostics

[www.operon.es](http://www.operon.es) - [sales@operon.es](mailto:sales@operon.es) - Tel.: +34 976 503 597 - Spain



Rev. 1

# Clostridium difficile

Name Nombre	Analyte Analito	Sample Muestra	Sensitivity Sensibilidad	Specificity Especificidad	Stick format	Simple format
GDH	Clostridium difficile	Stool samples Heces	96%	98%	X	X
2 A-B diff	Clostridium difficile	Stool samples Heces	94,9%	96,1%	X	X

## Sample preparation / Preparación de muestra



### Stick GDH



### Simple GDH



### Stick 2 A-B diff

**INDEPENDENT bands for toxin A and toxin B**



### Simple 2 A-B diff



# Infectious diseases *Enfermedades infecciosas*

## Norovirus, Rotavirus, Rota-Adeno, Rota-Noro

Immunochromatographic test for the diagnosis of Infectious diseases in stool samples  
*Test inmunocromatográfico para el diagnóstico de enfermedades infecciosas en heces*



[www.operon.es](http://www.operon.es) - [sales@operon.es](mailto:sales@operon.es) - Tel.: +34 976 503 597 - Spain



Rev. 4

## ROTAVIRUS, ROTA-ADENO

Name Nombre	Analyte Analito	Sample Muestra	Sensitivity Sensibilidad	Specificity Especificidad	Stick format	Simple format
<small>Reference method: ELISA</small>						
Rotavirus	Rotavirus	Stool samples	99.9%	98.6%	X	X
Rota-Adeno		Heeces	98.8%	92.2%		
Rota-Adeno	Adenovirus	Stool samples	97.6%	96.9%	X	X
		Heeces				


### SIMPLE ROTA, ROTA-ADENO




#### Results / Resultados: Rota - 5 min



#### Results / Resultados: Rota-Adeno - 10 min



### STICK ROTA, ROTA ADENO




#### Results / Resultados: Rota - 5 min



#### Results / Resultados: Rota-Adeno - 10 min



# NOROVIRUS



Name Nombre	Analyte Analito	Sample Muestra	Sensitivity Sensibilidad	Specificity Especificidad	Stick format	Simple format
Norovirus	Norovirus	Stool samples Heces	85%	97%	X	X

## Sample preparation / Preparación de muestra



## SIMPLE NOROVIRUS



## STICK NOROVIRUS

### INDEPENDENT bands for Norovirus GI and GII





# ROTA-NORO



Name Nombre	Analyte Analito	Sample Muestra	Sensitivity Sensibilidad	Specificity Especificidad	Stick format	Simple format
Rota-Noro	Rotavirus	Stool samples Heces	95%	97%	X	X
	Norovirus		85%	97%		

## Sample preparation / Preparación de muestra



## SIMPLE ROTA-NORO



## STICK ROTA-NORO

### INDEPENDENT bands for Norovirus and Rotavirus

